



# Clorox Healthcare Equipment Compatibility Guide



Healthcare facilities strive to reduce the risk of healthcare-associated infections (HAIs) associated with medical equipment. Yet at the same time, those devices must be safeguarded and maintained in good condition.

Medical equipment and devices are made of many different materials all of which may react differently to contact with disinfecting chemicals. Most equipment manufacturers supply their customers with cleaning & care guides that include information about what chemicals or products can be safely used on the equipment. Failure to follow these guides can result in voided warranties and expensive equipment replacements.

Clorox Healthcare® works directly with equipment manufacturers to test our products for compatibility with common medical equipment. Below is a list of equipment that includes approval to use Clorox Healthcare® products in the cleaning & care guide. If you are interested in using Clorox Healthcare® products on a piece of equipment that is not listed below, please contact your Clorox Healthcare® sales rep or email: [PPDTechnicalService@clorox.com](mailto:PPDTechnicalService@clorox.com).



For more information contact your Clorox sales representative or call **1.866.789.4973**  
Visit: [www.cloroxhealthcare.ca](http://www.cloroxhealthcare.ca)  
Email: [cleaning@clorox.com](mailto:cleaning@clorox.com)

© 2017 Clorox Professional Products Company, 150 Biscayne Crescent, Brampton, Ontario L6W 4V3

# Table of Contents

- 3 Beds/Pads/Vinyl
- 4 Blood Glucose Meters
- 5 Furniture/Upholstery/Apparel
- 7 Hemodialysis Systems
- 8 Infusion Pumps
- 9 Keyboards/Phones/Cases
- 10 Lights
- 11 Monitors/Consoles/Display/Transmitters
- 12 Radiology
- 13 Ultrasounds/Transducers/Probes
- 26 Ventilator Systems
- 27 Miscellaneous

## Endorsement Symbols



Product is mentioned by name in Cleaning & Care guide and is approved for use



Active ingredient/level is listed in Cleaning & Care guide as approved for use, but product is not listed by name









## Beds/Pads/Vinyl

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 <p>Clorox Healthcare® Professional Disinfecting Bleach Wipes</p>	 <p>Clorox Healthcare® Fuzion™ Cleaner Disinfectant</p>	 <p>Clorox® Disinfecting Bleach</p>	 <p>Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes</p>
<b>INFAB</b>				
X-ray aprons				<input checked="" type="checkbox"/>
<b>MIDMARK</b>				
Medical Exam Bed (Acute & Physicians)	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>
<b>MIZUHOSI</b>				
Surgical Pads				
<b>SKYTRON</b>				
Surgical Pads	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>

## Blood Glucose Meters

	Bleach			Hydrogen Peroxide
Manufacturer & Product	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>ABBOTT</b>				
FreeStyle InsuLinx®				
FreeStyle Lite®, Freedom Lite®			✓	
i-STAT® System			✓	
Precision Xtra™ Glucose & Ketone			✓	
<b>ARKRAY</b>				
ASSURE® 4			✓	
ASSURE® Platinum Glucose Meter	✓		✓	✓
ASSURE® Pro			✓	
<b>BAYER</b>				
CONTOUR® Link	✓			
CONTOUR® USB Meters			✓	
<b>EPS BIOTECH</b>				
EasyMax® Glucose Meter				✓
<b>LIFESCAN</b>				
OneTouch® Blood Glucose Meters			✓	

## Blood Glucose Meters (continued from previous page)

Manufacturer & Product	Bleach			Hydrogen Peroxide
				
	Clorox Healthcare® Professional Disinfecting Bleach Wipes	Clorox Healthcare® Fuzion™ Cleaner Disinfectant	Clorox® Disinfecting Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>NOVA BIOMEDICAL</b>				
StatStrip®				
<b>ROCHE</b>				
ACCU-CHEK® Advantage Meter System				
ACCU-CHEK® Inform II				
<b>TAIDOC</b>				
TaiDoc				
<b>MEDLINE</b>				
EVENCARE G2				
EVENCARE G3				

## Furniture/Upholstery/Apparel












Manufacturer & Product	Bleach			Hydrogen Peroxide
				
	Clorox Healthcare® Professional Disinfecting Bleach Wipes	Clorox Healthcare® Fuzion™ Cleaner Disinfectant	Clorox® Disinfecting Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>KWALU</b>				
All furniture products				
<b>HEALTHCENTRIC</b>				
IC+ Upholstery Solution™				
<b>PRIME MEDICAL</b>				
SAF-T™ Scrubs <sup>1</sup>				

<sup>1</sup> Clorox Healthcare® and Clorox® disinfectants are intended for use on hard, non-porous surfaces. SAF-T™ Scrubs endorsement covers color integrity under repeated contact with bleach disinfectants (e.g. user wearing scrubs while disinfecting surfaces)

## Hemodialysis Systems





Manufacturer & Product	Bleach			Hydrogen Peroxide
		 <p>Clorox Healthcare® Professional Disinfecting Bleach Wipes</p>	 <p>Clorox Healthcare® Fuzion™ Cleaner Disinfectant</p>	 <p>Clorox® Disinfecting Bleach</p>
Fresenius Medical Care				
2008K series			✓	
2008T series			✓	
Granuflo Concentrate Unit 1			✓	

# Infusion Pumps

Manufacturer & Product	Bleach			Hydrogen Peroxide
				
<b>BD Alaris (Carefusion)</b>				
CareFusion Alaris® Systems				
<b>Hospira</b>				
Omni-Flow 4000® Plus IV Medication Mgmt. System				
GemStar® Infusion Pumps				
<b>Agilia Intuitive Gen.</b>				
Volumetric Infusion Pump				
<b>Baxter</b>				
Sigma Spectrum				





## Keyboards/Phones/Cases

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>Whitney Medical</b>				
eShield® Screen Protectors and Cases	✓		✓	
<b>KeySource Int'l</b>				
San-a-Key® Keyboards	✓			✓
<b>Man and Machine</b>				
Keyboards and Mice	✓		✓	✓
<b>Otterbox</b>				
Phone and iPad® cases	✓			✓
<b>Seal Shield</b>				
Silver Stream™ Medical Grade Keyboards	✓		✓	✓








## Lights

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 <p>Clorox Healthcare® Professional Disinfecting Bleach Wipes</p>	 <p>Clorox Healthcare® Fuzion™ Cleaner Disinfectant</p>	 <p>Clorox® Disinfecting Bleach</p>	 <p>Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes</p>
<b>SKYTRON</b>				
<b>Surgical Lighting</b>				
<b>WELCH ALLYN</b>				
<b>Exam Lighting</b>				
<b>Green Series™ 300 General Exam Lights</b>				
<b>TRUMPF</b>				
<b>TruLight® iLED® Medical Lighting</b>				

## Monitors/Consoles/Display/Transmitters

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>BARCO</b>				
Eonis® Diagnostic Displays				
<b>GENERAL ELECTRIC COMPANY</b>				
Transport Pro Monitors				
<b>MIDMARK</b>				
IQvitals Monitors				
<b>WELCH ALLYN</b>				
Connex Vital Signs Monitor 600 Series				
<b>PHILLIPS</b>				
IntelliVue Patient Mx400 / 450 / 500 / 600 / 700 / 800				
IntelliVue Telemetry Transmitters				
IntelliVue Monitors				
Avalon Fetal Monitors				

## Radiology

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 <p>Clorox Healthcare® Professional Disinfecting Bleach Wipes</p>	 <p>Clorox Healthcare® Fuzion™ Cleaner Disinfectant</p>	 <p>Clorox® Disinfecting Bleach</p>	 <p>Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes</p>
<b>AGFA HEALTHCARE</b>				
<b>Computed Radiology Cassettes</b>				
<b>DR Detectors (Dx-D 10C, 10G, 20C, 20G)</b>				
<b>HOLOGIC</b>				
<b>Mammography</b>				

## Ultrasounds/Transducers/Probes

Manufacturer & Product	Bleach			Hydrogen Peroxide
		 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach
GENERAL ELECTRIC COMPANY				
LOGIQ™ Transducers/ Ultrasounds				
10L, 10LB, 10LB-RS, 10S, 10S-RS,				
11L, 11L-D				
12L, 12L-RS, 12L-SC				
125-D				
16-RS				
3.5C, 3.5CS, 3.5 gray lens, 3C, 3Cb, 3CRF, 3CRF-D, 3C-RS, 3S, 3SC-RS, 3S-D, 3S- LC, 3SR, 3S-RC, 3S-RS, 3S-SC, 3V, 3V-D				
4C				
4C-A, 4C-D				
4C-RS, 4C-SC				

## Ultrasounds/Transducers/Probes (continued from previous page)

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>GENERAL ELECTRIC COMPANY cont'd</b>				
4D10L				
4D16L, 4D3C-L, 4D5C-L, 4D8C, 4DE7C, 4S, 4V-D	✔			
546L, 5C, 5CS, 5S, 5S-LC,	✔			
5S-RS	✔			
6S, 6S-D	✔			
739L, 7L, 7S, 7S-R5	✔			
8C, 8C-RS	✔			
9L, 9L-D, 9L-RS	✔			
AB2-5, AB2-7, AB2-7-D, AB2-7-RS, AB4-8	✔			
AC2-5	✔			
BE9C, BE9CS, BE9CS-D, BE9CS-RS	✔			

## Ultrasounds/Transducers/Probes (continued from previous page)











	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>GENERAL ELECTRIC COMPANY cont'd</b>				
C1-5-RS, C1-6-D C1-6VN-D, C31, C3-10-D, C358, C36, C364, C386, C4-8-D, C55, C551, C721, CAE, CBF, CS, CZB				
C2-5-RC, C2-9-D				
Dual Probe Vscan				
E72, E721, E8C, E8C-RS, E8CS, E8CS-RC, E8CS-RS, E8CS-SC, ERB				
eM6C				
G3C Probe on Vscan, C3S Probe on Vscan				
i12L, i12L-RS				
I739, I739-RS, IC5-9, IC5-9-D, IC5-9H, IC5-9W-RS, IC9-RS				
L12n-SC, L3-12D, L39, L4-12t-RS, L6-12-RS, L76, L764, L8-18i-D, L8- 18i-RS, L8-18i-SC, LA39, LB, LD, LE, LH, LI, LT				

## Ultrasounds/Transducers/Probes (continued from previous page)

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>GENERAL ELECTRIC COMPANY cont'd</b>				
M12L, M12L-H, M3S, M4S, M4S-D, M5Sc-RS, M6C, M7C, M7C-H, ML6-15, ML6-15-D, ML6-15-RS, MTZ, MZ				
M3-8				
PA2-5P, PCW4.0				
RAB2-5, RAB2-5-D, RAB2-5L, RAB2-5-RS, RAB4-8-D, RAB4-8L, RAB4-8P, RAB4-8-RS, RAB6-D, RAB6-RS, RAM3-8, RIC5-9, RIC5-9A-RS, RIC5-9-D, RIC5-9H, RIC5-9-RS, RIC5-9W, RIC5-9W-RS, RIC6-12-D, RM14L, RM6C, RNA5-9, RNA5-9-D, RNA5-9-RS, RRE5-10-D, RRE6-10, RRE6-10-D, RSM5-14, RSP5-12, RSP6-12, RSP6-16, RSP6-16-D, RSP6-16-RS				
S220, S222, S316, S317, S611, SCW2.0,				
SP10-16-RS, SP4-10, SP6-12				
T739, T739-RS				



## Ultrasounds/Transducers/Probes (continued from previous page)

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>PHILLIPS</b>				
<b>Ultrasound Transducers &amp; Probes (SystemNames):</b> <b>T = Transducer / C = Cable and Connector / H = Handle and Steering Mechanism</b>				
3D6-2 (HD11), 3D8-4, 3D9-3v (HD11, HD15, iU22) ,		 T,C		 T,C
7.5MHZ Endo		 T,C		 T,C
C5-2 (HD5, HD6, HD7, HD11, HD15, iE33, iU22) , C6-2 Compact (CX30, Sparq), C6-3 , C8-4v, C8-5 , C8-5 Compact (CX30, CX50) , C9-3io , C9-4, C9-4ec, C9-5ec		 T,C		 T,C
L8-4, L9-3 , L11-3, L12-3 (HD5, HD6, HD7, HD11, HD15), L12-4 Compact (CX30), L12-5 50 Compact (CX50), L12-5 50 mm, L15-7io (HD11, HD15, iE33, iU22), L17-5		 T,C		 T,C

## Ultrasounds/Transducers/Probes (continued from previous page)

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>PHILLIPS cont'd</b>				
<b>Ultrasound Transducers &amp; Probes (SystemNames):</b> <b>T = Transducer / C = Cable and Connector / H = Handle and Steering Mechanism</b>				
S3-1, S4-1 (iU22), S4-2, S4-2 Compact (CX30, Sparq), S5-1 (CX50, iE33, iU22), S5-2, S8, S8-3, S8-3 Compact (CX30, CX50, HD5), S12, S12-4, S12-4 Compact (CX50)		 T,C		 T,C
V6-2 (HD11XE, HD15, iU22), V8-4		 T,C		 T,C
X3-1, X5-1, X7-2 (iE33, iU22)		 T,C		 T,C
3D6-2 (iU22), 3D9-3v (Affiniti, ClearVue, EPIQ)		 T,C		 T,C
BP10-5ec		 T,C		 T,C
C5-1 (Affiniti, EPIQ), C5-2 - ClearVue, C6-2 Compact (Affiniti), C8-4v, C8-5 Compact (Affiniti, EPIQ), C9-2, C9-3v, C9-4v, C9-4v Compact, C10-3v, C10-4ec		 T,C		 T,C





## Ultrasounds/Transducers/Probes (continued from previous page)

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>PHILLIPS cont'd</b>				
<b>Ultrasound Transducers &amp; Probes (SystemNames):</b> <b>T = Transducer / C = Cable and Connector / H = Handle and Steering Mechanism</b>				
L12-3 (Affiniti, CX50, EPIQ), L12-4 X7-2 (EPIQ), L12-4 Compact (Affiniti), L12-5 50 Compact (Affiniti, EPIQ), L15-7io (Affiniti, CX30, CX50, EPIQ), L15-7io (HD11, HD15, iE33, iU22), L18-5		 T,C		 T,C
S4-1 (ClearVue), S4-2 Compact (Affiniti), S5-1 (Affiniti, EPIQ), S8-3 Compact (Affiniti, EPIQ), S12-4 Compact (Affiniti, EPIQ)		 T,C		 T,C
V6-2 (Affiniti, ClearVue, EPIQ, HD5), VL13-5, VL13-5 Compact		 T,C		 T,C
X5-1 (EPIQ), X6-1		 T,C		 T,C
1.9MHz CW Pencil, 1.9MHz TCD		 T,C		 T,C
5.0MHz CW Pencil		 T,C		 T,C

## Ultrasounds/Transducers/Probes (continued from previous page)

Manufacturer & Product	Bleach			Hydrogen Peroxide
				
	Clorox Healthcare® Professional Disinfecting Bleach Wipes	Clorox Healthcare® Fuzion™ Cleaner Disinfectant	Clorox® Disinfecting Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>PHILLIPS cont'd</b>				
<b>Ultrasound Transducers &amp; Probes (SystemNames):</b>				
<b>T = Transducer / C = Cable and Connector / H = Handle and Steering Mechanism</b>				
D2cwc, D2tcd, D5cwc		 T,C		 T,C
TCD		 T,C		 T,C
S6-2mpt S7-3t, S7-2t Omni, S8-3t		 C		 T,C,H
OMNI III TEE		 C		 T,C,H
X7-2t		 C		 T,C,H
L10-4lap		 T,C		 T,C,H
15-6L		 T,C		 C

## Ultrasounds/Transducers/Probes (continued from previous page)

Manufacturer & Product	Bleach			Hydrogen Peroxide
				
	Clorox Healthcare® Professional Disinfecting Bleach Wipes	Clorox Healthcare® Fuzion™ Cleaner Disinfectant	Clorox® Disinfecting Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>PHILLIPS cont'd</b>				
<b>CO2 Transducer</b>				
M1460A			✓	
Avalon Fetal Transducer			✓	
<b>SIEMENS</b>				
Acuson CV70 , Acuson SC2000 , Acuson Cypress , Acuson X500 , Acuson X300, Acuson Antares , Acuson X150, and Sonoline G20 UltraSound Systems				
3V2c, 4V1c, 5.0L45, 5.0P10, 6C1HD, 7.5L75S, 8V3, 10V4	✓			
BE9-4	✓			
C6-2, C7F2, C8-5, CW2, CW5, C5-2, C5F1, CH4-1, CH6-2, CX5-2	✓			
EC9-4, EV9-4	✓			
L10-5	✓			

## Ultrasounds/Transducers/Probes (continued from previous page)

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>SIEMENS cont'd</b>				
MPT7-4	✓			
P4-2, P9-4, P10-4, PH4-1, PX4-1	✓			
AUX CW	✓			
VF13-5-SP, VF7-3, VFX13-5, V5Ms, VF13-5-SP, V7M	✓			
<b>FUJIFILM SONOSITE, INC.</b>				
<b>System and Mini-Dock</b>				
Edge™				✓
NanoMaxx®				✓
S Series™				✓
M-Turbo® M-OB/GYN Office™				✓

## Ultrasounds/Transducers/Probes (continued from previous page)





Manufacturer & Product	Bleach			Hydrogen Peroxide
				
	Clorox Healthcare® Professional Disinfecting Bleach Wipes	Clorox Healthcare® Fuzion™ Cleaner Disinfectant	Clorox® Disinfecting Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>FUJIFILM SONOSITE, INC. cont'd</b>				
MicroMaxx®				
TITAN®				
Sonosite® 180PLUS™				
iLook®15, iLook®25				
<b>Nanosonics Trophon EPR</b>				
C11, C11e, C11x, C11 (Cast Urethane), C15 (Cast Urethane)				
C11n				
C15, C15e				
C60				
C60e, C60x, C60n				
D2, D2x				

## Ultrasounds/Transducers/Probes (continued from previous page)

Manufacturer & Product	Bleach			Hydrogen Peroxide
				
	Clorox Healthcare® Professional Disinfecting Bleach Wipes	Clorox Healthcare® Fuzion™ Cleaner Disinfectant	Clorox® Disinfecting Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
FUJIFILM SONOSITE, INC. cont'd				
HFL38, HFL38x				
HFL50x				
L25, L25e, L25x				
L25n				
L38, L38e, L38x, L52, L52e (Vet Only)				
L38n, L38xi, L52x (Vet Only), L52n (Vet Only)				
P10, P10x, P17, P21x				
P21n				
C8, C8e, C8x				
ICT, ICTx				





## Ultrasounds/Transducers/Probes (continued from previous page)

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 <p>Clorox Healthcare® Professional Disinfecting Bleach Wipes</p>	 <p>Clorox Healthcare® Fuzion™ Cleaner Disinfectant</p>	 <p>Clorox® Disinfecting Bleach</p>	 <p>Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes</p>
<b>FUJIFILM SONOSITE, INC. cont'd</b>				
ICTn	<input checked="" type="checkbox"/>			
SLA, SLAx, SLT	<input checked="" type="checkbox"/>			
<b>Accessories</b>				
NanoMaxx® Dock	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>
PowerPack	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>
Power Supply P06936, P08792, P09823, P09714	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>
<b>Welch Allyn</b>				
SureTemp® Therms and Probes	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	





## Ventilator Systems

Manufacturer & Product	Bleach			Hydrogen Peroxide
		 <p>Clorox Healthcare® Professional Disinfecting Bleach Wipes</p>	 <p>Clorox Healthcare® Fuzion™ Cleaner Disinfectant</p>	 <p>Clorox® Disinfecting Bleach</p>
<b>Covidien</b>				
<b>Puritan Bennett</b>				
540 series (409 was OK'd for touchscreen and arm)			✓	
840 series (409 was OK'd for touchscreen and arm)			✓	
980 series (409 was OK'd for touchscreen and arm)			✓	
<b>Kendall SCD</b>				
700 Sequentia Compression System				
<b>AV Impulse</b>				
6060 Controller and Accessories				
<b>Phillips</b>				
Trilogy 100 System			✓	

## Miscellaneous

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>Blickman</b>				
Surgical Steel				
<b>Carefusion</b>				
Pyxis® Supply Systems				
<b>Enovate Medical</b>				
Computer Carts (EMC)				
<b>Midmark</b>				
Upholstery				
Manual Examination Tables				
IQvitals® Blood Pressure Cuffs				
<b>Roche</b>				
Cobas® Urisys 1100® Urine Analyzers				
<b>Sechrist</b>				
Chambers				

## Miscellaneous (continued from previous page)

	Bleach			Hydrogen Peroxide
<b>Manufacturer &amp; Product</b>	 Clorox Healthcare® Professional Disinfecting Bleach Wipes	 Clorox Healthcare® Fuzion™ Cleaner Disinfectant	 Clorox® Disinfecting Bleach	 Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
<b>Welch Allyn®</b>				
<b>Accessory Cable Management Mobile Stand 4800/4900 Series</b>				
<b>Propaq Encore® Blood Pressure Cuffs</b>				
<b>Harvey™ Elite® Stethoscopes</b>				
<b>LumiView™ Binocular Microscope</b>				
<b>Digital MacroView™ Otoscopes and Halogen HPX® Illuminators</b>				
<b>RF Surgical Systems</b>				
<b>RF Assure Detection System</b>				